




USER MANUAL for Build v0.1.4

TABLE OF CONTENTS

INTRODUCTION	1
OVERALL WORKFLOW	2
FUNCTIONAL MODULES.....	3
SOFTWARE AVAILABILITY.....	4
INSTALLATION INSTRUCTIONS.....	4
MODULE 1: PROTEIN IDENTIFICATION	5
MODULE 2: QUANTITATIVE FEATURE GENERATION.....	7
MODULE 3: DATA INTEGRATION AND QUANTIFICATION	11
MODULE 4: POST-QUANTIFICATION DATA PROCESSING.....	17
CONTACT INFORMATION.....	18

INTRODUCTION




IonStar is an MS1-based quantitative method for label-free proteomics experiments, devised to address issues related with quantitative precision, missing data, and false-positive discovery of protein changes in large-cohort analysis. IonStar comprises of two parts: experimental procedures (left panel) and a proteomics data analysis pipeline (right panel). Details of the experimental procedures can be found in the following literature:

- ◆ Shen X, Shen S, Li J, Hu Q, Nie L, Tu C, Wang X, Orsburn B, Wang J, Qu J*. An IonStar experimental strategy for MS1 ion current-based quantification using ultra-high-field Orbitrap: reproducible, in-depth and accurate protein measurement in larger cohorts. *J Proteome Res.* 16(7):2445-2456. (2017)
- ◆ An B, Zhang M, Johnson RW, Qu J*. Surfactant-Aided Precipitation/On-Pellet-Digestion (SOD) Procedure Provides Robust and Rapid Sample Preparation for Reproducible, Accurate and Sensitive LC-MS Quantification of Therapeutic Protein in Plasma and Tissues. *Anal Chem.* 87(7):4023-9 (2015)

This manual will focus on the data analysis pipeline part of IonStar, including the functional modules, overall workflow, and the detailed usage of each module. An empirical example will also be provided to help IonStar users to run the pipeline in their own environment.

OVERALL WORKFLOW



As shown in the figure above, LC-MS raw files are processed in two separated routes to generate protein identification and quantification results. In one route, rawfiles are imported into database search engines for peptide identification, and generated search results are processed for FDR control, peptide inference, and protein group. The final output is a **.csv spectrum report** containing the filtered PSM list. In the other route, SIEVETM is used for quantitative feature generation, including ChromAlign-based RT correction and direct ion-current extraction (DICE). The features are stored in an **.sdb SQL database**. The spectrum report and feature database are merged by a customized R package, IonStarStat, to match features with corresponding PSMs. The IonStarStat packages also serves to perform dataset-wide normalization, multivariate algorithm-based outlier rejection, and data aggregation.


FUNCTIONAL MODULES

Module 1: Protein Identification


General workflow




MS-GF+-based workflow



Module 2: Quantitative feature extraction



Module 3: Data integration and quantification



Module 4: Post-quantification data processing



Processing of quantitative results can be achieved by StarGazer, an R-Shiny App currently under development. StarGazer will be available in the next build of IonStar.

SOFTWARE AVAILABILITY

SIEVE™ v2.2 SP2 is a commercial software from Thermo Fisher Scientific. Please contact Thermo Fisher Scientific regarding the quote for SIEVE™. To ensure of proper performance of SIEVE™, we recommend running SIEVE™ on a PC with at least 16-core processors and at least 192 GB RAM.

R package IonStarStat and related scripts (IonStarSPG.R, Load.R, Run IonStar code.R) can be downloaded from <https://github.com/shichens1989/IonStarStat>. All operations in this manual are accomplished under R version 3.2.4 and Rstudio ver1.0.136.

INSTALLATION INSTRUCTIONS

IonStarStat package can be installed directly in RStudio by running the following command in the R Console:

```
#####Install dependency and ion star(the beta version called MStest)#####  
#Install dependencies "RSQLite""MCMCglmm""affyPLM""mvoutliers"  
source("https://bioconductor.org/biocLite.R")  
biocLite("affyPLM")  
biocLite("MCMCglmm")  
biocLite("RSQLite")  
biocLite("mvoutliers")  
Install.packages(("File Directory\\IonStarStat_0.1.4.tar.gz"))
```

Alternatively, the **Install Packages** under the **Tools** menu can be used for installation. Yet the supporting packages need to be installed manually as well.

MODULE 1: PROTEIN IDENTIFICATION

The first input into the IonStarStat package is a spectrum report containing a filtered PSM list. This report can be exported conveniently from a number of softwares, e.g. Proteome Discoverer, Scaffold. Here we will demonstrate the use of the IonStarSPG.R spectrum report generator using results obtained from MS-GF+ (database search) and IDPicker (post-search processing).


I. Input files

Search results .tsv files (**Sample name.raw.tsv**) converted from .mzid files using MSGFPlus.jar (-showQValue = 1, -shoeDecoy = 1, -unroll = 1)

Protein list .tsv file exported from the Protein View tab of the .idpDB file

Peptide list .tsv file exported from the Peptide View tab of the .idpDB file

Spectrum list .tsv file exported from the Spectrum View tab of the .idpDB file (Tree grouping = Spectrum)



II. Spectrum report generation with IonStarSPG.R

1. Confirm that all input files are located under the same directory.

2. Generate spectrum report

a. Change the input and output file names:

Line 4~7

```
4 setwd("File Directory")
5 protein <- read.csv("Protein list.tsv",header=TRUE,sep="\t",fileEncoding="windows-1252")
6 peptide <- read.table("Peptide list.tsv",header=TRUE,sep="\t",fileEncoding="windows-1252")
7 spectrum <- read.table("Spectra list.tsv",header=TRUE,sep="\t",fileEncoding="windows-1252")
```

Line 70

```
70 write.csv(spectrumrepo_total, "Spectrum report.csv", row.names=FALSE)
```

b. Run the script by **Source** function.

c. The output (Spectrum report.csv) contains the filtered PSM from each LC-MS run as follows:

#	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA				
1	SpecFile	SpecID	ScanNum	FragMeth	Precursor	IsotopeEn	Precursor	Charge	DeNovoSc	MSGFScor	SpecEVal	LEValue	QValue	PepQVal	Peptide	S	Accession	Cluster	Count	Coverage	Protein	GI	Distinct	Pi	Distinct	N	Filtered	S	Descriptive	Accession	2
2	171103_IAScan=2719	27198	HCD		612.5788	1	2.617393	4	166	163	9.65E-29	1.80E-21	0	0	IVVVVVVG	sp	Q61201	1	1	58.18966	1	8	8	136	Platelet-a	sp	Q61205	PA1B3_MOUSE			
3	171103_IAScan=1134	11340	HCD		1134.528	1	-3.84731	2	165	162	1.34E-25	2.51E-18	0	0	SGEGENP	sp	Q61201	1	1	58.18966	1	8	8	136	Platelet-a	sp	Q61205	PA1B3_MOUSE			
4	171103_IAScan=5277	52776	HCD		483.329	0	0.126281	2	111	111	3.60E-11	0.000656	0.000104	0.00015	VVVVGLP	sp	Q61201	1	1	58.18966	1	8	8	136	Platelet-a	sp	Q61205	PA1B3_MOUSE			
5	171103_IAScan=1136	11364	HCD		1134.029	0	-1.07643	2	189	175	8.59E-25	1.60E-17	0	0	SGEGENP	sp	Q61201	1	1	58.18966	1	8	8	136	Platelet-a	sp	Q61205	PA1B3_MOUSE			
6	171103_IAScan=4358	43586	HCD		867.5212	1	-0.10438	2	138	99	2.79E-12	5.19E-05	0	0	LPHLTLTSS	sp	Q8BK01	2	1	16.31579	2	2	2	16	Transmem	sp	Q8BK08	TMM11_MOUSE			
7	171103_IAScan=5969	59690	HCD		893.419	0	-12.9624	2	109	61	4.40E-12	8.14E-05	0	0	EFQEGPLL	sp	Q8C5L1	3	1	10.47009	3	3	3	58	Inositol	sp	Q8C5L6	INPSK_MOUSE			
8	171103_IAScan=5693	56931	HCD		861.4202	0	-1.06281	2	149	125	2.04E-15	3.79E-08	0	0	IQCSIMGLT	sp	Q08791	4	1	4.41989	4	2	2	15	Eukaryoti	sp	Q08796	EF2K_MOUSE			
9	171103_IAScan=1284	12843	HCD		616.3039	0	0.396137	2	115	109	5.20E-14	9.55E-07	0	0	AGGFTLGT	sp	P11672	6	1	25	6	5	5	132	Neutroph	sp	P11672	NGAL_MOUSE			
10	171103_IAScan=4309	43094	HCD		702.8854	0	0	2	139	134	3.41E-18	6.39E-11	0	0	WVVVGLA	sp	P11672	6	1	25	6	5	5	132	Neutroph	sp	P11672	NGAL_MOUSE			
11	171103_IAScan=2348	23484	HCD		431.2393	-1	-4.03146	2	69	69	1.12E-10	0.002022	0.000277	0.000401	IADFGGLAR	sp	P39686	8	1	21.04283	8	9	9	141	Tyrosine	sp	P39688	FYN_MOUSE			
12	171103_IAScan=7987	79877	HCD		612.3014	0	4.286326	2	115	115	2.09E-13	9.83E-06	0	0	UEDENETV	sp	P39686	8	1	21.04283	8	9	9	141	Tyrosine	sp	P39688	FYN_MOUSE			
13	171103_IAScan=3083	30834	HCD		618.3144	1	-2.71508	2	122	120	2.84E-15	5.22E-08	0	0	WTAPEAR	sp	P39686	8	1	21.04283	8	9	9	141	Tyrosine	sp	P39688	FYN_MOUSE			
14	171103_IAScan=2488	24883	HCD		431.7436	0	2.120539	2	103	103	2.23E-10	0.004044	0.000592	0.000401	IADFGGLAR	sp	P39686	8	1	21.04283	8	9	9	141	Tyrosine	sp	P39688	FYN_MOUSE			
15	171103_IAScan=2388	23883	HCD		774.0774	0	1.103886	3	179	105	9.05E-14	1.69E-06	0	0	MERPGRG	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
16	171103_IAScan=3101	31011	HCD		1054.183	0	1.042167	3	243	93	1.40E-13	2.63E-06	0	0	SSVVDMM	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
17	171103_IAScan=5282	52828	HCD		888.4451	0	0.961785	3	244	107	5.79E-13	1.08E-05	0	0	AQVLDCC	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
18	171103_IAScan=3096	30962	HCD		632.9118	0	-0.09644	5	215	162	2.46E-22	4.61E-15	0	0	SSVVDMM	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
19	171103_IAScan=5066	50669	HCD		1388.681	0	2.461313	2	189	162	8.34E-22	1.56E-14	0	0	TMSIVSVY	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
20	171103_IAScan=8006	80069	HCD		674.3995	0	2.53409	2	131	126	7.32E-15	1.35E-07	0	0	SVLVDFLI	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
21	171103_IAScan=5283	52834	HCD		666.5853	0	0.549384	4	235	206	4.48E-22	8.36E-15	0	0	AQVLDCC	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
22	171103_IAScan=5302	53028	HCD		666.3282	-1	-8.90538	4	200	130	2.17E-17	4.04E-10	0	0	AQVLDCC	sp	Q9JHU1	9	1	42.90844	9	17	17	549	Inositol-3	sp	Q9JHU9	INO1_MOUSE			
23	171103_IAScan=4952	49525	HCD		853.4715	0	0.786655	2	188	183	6.51E-19	1.21E-11	0	0	NLNQALLI	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
24	171103_IAScan=6369	63698	HCD		820.899	1	1.506723	4	264	261	2.12E-36	9.97E-29	0	0	ALFQDVQ	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
25	171103_IAScan=4946	49467	HCD		569.3164	0	0.214416	3	134	130	2.82E-19	5.24E-12	0	0	NLNQALLI	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
26	171103_IAScan=2634	26348	HCD		938.467	0	1.690967	2	156	150	1.77E-19	3.29E-12	0	0	ALFQDVQ	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
27	171103_IAScan=5168	51685	HCD		791.7334	0	1.002178	3	179	165	4.42E-22	8.25E-15	0	0	VAGPQPA	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
28	171103_IAScan=2887	28874	HCD		711.8348	0	1.45764	2	135	132	4.00E-17	7.38E-10	0	0	TQEAMEA	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
29	171103_IAScan=5162	51628	HCD		1187.596	1	-0.59039	2	250	236	2.68E-24	5.01E-17	0	0	VAGPQPA	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
30	171103_IAScan=6375	63755	HCD		1094.193	1	-1.6547	3	255	212	1.70E-25	3.19E-18	0	0	ALFQDVQ	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
31	171103_IAScan=6314	63148	HCD		528.2744	0	-0.11554	2	98	97	5.05E-13	9.21E-06	0	0	MGNHLT	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
32	171103_IAScan=1921	19214	HCD		740.8408	0	-19.1133	2	117	75	1.72E-10	0.003182	0.000422	0.000114	QNVSTVE	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
33	171103_IAScan=1856	18569	HCD		494.2394	0	0.740959	3	106	76	2.91E-11	0.000537	7.90E-05	0.000114	QNVSTVE	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
34	171103_IAScan=5217	52172	HCD		1187.596	1	-0.59039	2	250	243	6.38E-25	1.19E-17	0	0	VAGPQPA	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
35	171103_IAScan=4947	49476	HCD		569.3166	0	0.536039	3	167	161	6.39E-19	1.19E-11	0	0	NLNQALLI	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
36	171103_IAScan=6363	63637	HCD		1093.855	0	-4.91022	3	128	93	5.09E-12	9.53E-05	0	0	ALFQDVQ	sp	P29391	10	1	65.02732	10	13	13	547	Ferritin	sp	P29391	FRIL1_MOUSE			
37	171103_IAScan=2634	26348	HCD		938.467	0	1.690967	2	156	150	1.77E-19	3.29E-12	0	0	ALFQDVQ	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
38	171103_IAScan=1856	18569	HCD		494.2394	0	0.740959	3	106	76	2.91E-11	0.000537	7.90E-05	0.000114	QNVSTVE	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
39	171103_IAScan=4946	49467	HCD		569.3164	0	0.214416	3	134	130	2.82E-19	5.24E-12	0	0	NLNQALLI	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
40	171103_IAScan=4947	49476	HCD		569.3166	0	0.536039	3	167	161	6.39E-19	1.19E-11	0	0	NLNQALLI	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
41	171103_IAScan=4952	49525	HCD		853.4715	0	0.786655	2	188	183	6.51E-19	1.21E-11	0	0	NLNQALLI	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
42	171103_IAScan=6314	63148	HCD		528.2744	0	-0.11554	2	98	97	5.05E-13	9.21E-06	0	0	MGNHLT	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
43	171103_IAScan=1921	19214	HCD		740.8408	0	-19.1133	2	117	75	1.72E-10	0.003182	0.000422	0.000114	QNVSTVE	sp	P49946	10	1	30.60109	11	7	7	266	Ferritin	sp	P49946	FRIL2_MOUSE			
44	171103_IAScan=2381	23812	HCD		687.3893	0	2.752582	2	148	131	4.41E-13	8.13E-06	0	0	IEVEKPF	sp	P24363	11	1	63.88889	12	20	20	979	Peptidyl	sp	P24363	PP1B_MOUSE			
45	171103_IAScan=4638	46383	HCD		729.8593	1	-13.0948	2	155	130	3.55E-14	6.56E-07	0	0	DTNGSQFI	sp	P24363	11	1	63.88889	12	20	20	979	Peptidyl	sp	P24363	PP1B_MOUSE			


MODULE 2: QUANTITATIVE FEATURE GENERATION

Quantitative feature generation in IonStar is accomplished by SIEVE™ v2.2 SP2 (Thermo Scientific), which integrates ChromAlign for global 3-D chromatographic alignment and a DICE method for feature extraction.


I. Setting up the method

1. Load rawfiles into SIEVE™



To start the quantitative feature generation analysis, open SIEVE™ and select **File → Create new experiment**. On the **Designate Experiment Type** page, select the Experiment Type based on the study. For a case-control experiment, use **Two Sample Differential Analysis**; for multi-condition experiment (≥ 3 conditions including control), use **Control Compare Trend**.





Drag all rawfiles onto the **Raw File Selection** page.




Assign conditions to each file imported. For **Two Sample Differential Analysis**, assign Condition A and B in the two boxes; for **Control Compare Trend**, assign all conditions in the upper box and the control condition in the lower box.



Setup the method parameters. The parameters that needs to be modified (for LC-MS rawfiles acquired on a Thermo Orbitrap instrument under 120K MS1 resolution) are shown as below. All other parameters follow the default settings.



After loading files with the given parameters, For IonStar, users do not need to run the Identify process. MaxThreads also needs to be modified according to the configuration of the computer used for SIEVE™. For example, 6~8 threads are recommended for a computer with 16-core processors and 92 GB RAM. A suggested sequence on the **Process** tab is highlighted in the following figure.



After the Frame process, an .sdb file will be generated containing all quantitative features. For more detailed information about the use of SIEVE™, please refer to SIEVE™ User Guide:

<https://tools.thermofisher.com/content/sfs/manuals/Man-XCALI-97696-SIEVE-22-User-ManXCALI97696-A-EN.pdf>

MODULE 3: DATA INTEGRATION AND QUANTIFICATION

After protein identification and quantitative feature generation, the customized R package IonStarStat plus associating R scripts will be utilized to generate the final quantitative results. Steps involved in this module include:

- ◆ Generation of the annotated frame list by integrating the spectrum report and the SIEVE database;
- ◆ Removal of redundant frames;
- ◆ Aggregation of frame-level data to peptide level & inter-sample normalization of quantitative intensities;
- ◆ Multivariate mean variation-based outlier peptide detection;
- ◆ Shared peptide removal (optional);
- ◆ Peptide-to-protein aggregation via two different algorithms.

The instructions for each step using IonStarStat will be elaborated as follows.

I. Generation of the annotated frame list

1. Required packages for frame list annotation

To generate the frame list, the following R packages need to be installed: **XLConnect**, **RSQLite**.

2. Specify the input and output files

Input files: Line 1~5

```
1 setwd("File directory")
2 ##File name of SIEVE file
3 db<-"SIEVE database.sdb"
4 ##File name of Spectrum report
5 xls<-"Spectrum report.csv"
```

Output files: Line 37~38 or Line 71~72

```
37 write.csv(sub_data,"Annotated Frame List.csv",row.names=FALSE)
38 write.csv(colnames(sub_data),"Sample List.csv")
```

3. Generate the annotated frame list & the sample list

Generate the annotated frame list using the script.

Line 7~37 summarize scripts used for spectrum reports exported from Scaffold.

Line 39~70 summarize scripts used for spectrum reports merged by IonStarSPG.R.

Note that Line 36 or 70 needs to be modified.

```
36 sub_data<-FS[, "columns"]
```

The columns to be included are **Peptide AC**, **Peptide sequence**, **Frame ID**, and **Frame intensity in each sample**.

- ◆ If other software combinations are used to generate the spectrum report, the scripts can be modified so that column headers are consistent between the spectrum report and SIEVE database for matching (recommended for users familiar with R language).

An example of the exported Frame list is shown as follows:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE	AF	AG	AH	AI	AJ	AK	AL	AM	AN	AO	AP	AQ	AR	AS	AT	AU	AV	AW	AX	AY	AZ																																																
1	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
2	101	102	103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136	137	138	139	140	141	142	143	144	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159	160	161	162	163	164	165	166	167	168	169	170	171	172	173	174	175	176	177	178	179	180	181	182	183	184	185	186	187	188	189	190	191	192	193	194	195	196	197	198	199	200
3	201	202	203	204	205	206	207	208	209	210	211	212	213	214	215	216	217	218	219	220	221	222	223	224	225	226	227	228	229	230	231	232	233	234	235	236	237	238	239	240	241	242	243	244	245	246	247	248	249	250	251	252	253	254	255	256	257	258	259	260	261	262	263	264	265	266	267	268	269	270	271	272	273	274	275	276	277	278	279	280	281	282	283	284	285	286	287	288	289	290	291	292	293	294	295	296	297	298	299	300
4	301	302	303	304	305	306	307	308	309	310	311	312	313	314	315	316	317	318	319	320	321	322	323	324	325	326	327	328	329	330	331	332	333	334	335	336	337	338	339	340	341	342	343	344	345	346	347	348	349	350	351	352	353	354	355	356	357	358	359	360	361	362	363	364	365	366	367	368	369	370	371	372	373	374	375	376	377	378	379	380	381	382	383	384	385	386	387	388	389	390	391	392	393	394	395	396	397	398	399	400
5	401	402	403	404	405	406	407	408	409	410	411	412	413	414	415	416	417	418	419	420	421	422	423	424	425	426	427	428	429	430	431	432	433	434	435	436	437	438	439	440	441	442	443	444	445	446	447	448	449	450	451	452	453	454	455	456	457	458	459	460	461	462	463	464	465	466	467	468	469	470	471	472	473	474	475	476	477	478	479	480	481	482	483	484	485	486	487	488	489	490	491	492	493	494	495	496	497	498	499	500
6	501	502	503	504	505	506	507	508	509	510	511	512	513	514	515	516	517	518	519	520	521	522	523	524	525	526	527	528	529	530	531	532	533	534	535	536	537	538	539	540	541	542	543	544	545	546	547	548	549	550	551	552	553	554	555	556	557	558	559	560	561	562	563	564	565	566	567	568	569	570	571	572	573	574	575	576	577	578	579	580	581	582	583	584	585	586	587	588	589	590	591	592	593	594	595	596	597	598	599	600
7	601	602	603	604	605	606	607	608	609	610	611	612	613	614	615	616	617	618	619	620	621	622	623	624	625	626	627	628	629	630	631	632	633	634	635	636	637	638	639	640	641	642	643	644	645	646	647	648	649	650	651	652	653	654	655	656	657	658	659	660	661	662	663	664	665	666	667	668	669	670	671	672	673	674	675	676	677	678	679	680	681	682	683	684	685	686	687	688	689	690	691	692	693	694	695	696	697	698	699	700
8	701	702	703	704	705	706	707	708	709	710	711	712	713	714	715	716	717	718	719	720	721	722	723	724	725	726	727	728	729	730	731	732	733	734	735	736	737	738	739	740	741	742	743	744	745	746	747	748	749	750	751	752	753	754	755	756	757	758	759	760	761	762	763	764	765	766	767	768	769	770	771	772	773	774	775	776	777	778	779	780	781	782	783	784	785	786	787	788	789	790	791	792	793	794	795	796	797	798	799	800
9	801	802	803	804	805	806	807	808	809	810	811	812	813	814	815	816	817	818	819	820	821	822	823	824	825	826	827	828	829	830	831	832	833	834	835	836	837	838	839	840	841	842	843	844	845	846	847	848	849	850	851	852	853	854	855	856	857	858	859	860	861	862	863	864	865	866	867	868	869	870	871	872	873	874	875	876	877	878	879	880	881	882	883	884	885	886	887	888	889	890	891	892	893	894	895	896	897	898	899	900
10	901	902	903	904	905	906	907	908	909	910	911	912	913	914	915	916	917	918	919	920	921	922	923	924	925	926	927	928	929	930	931	932	933	934	935	936	937	938	939	940	941	942	943	944	945	946	947	948	949	950	951	952	953	954	955	956	957	958	959	960	961	962	963	964	965	966	967	968	969	970	971	972	973	974	975	976	977	978	979	980	981	982	983	984	985	986	987	988	989	990	991	992	993	994	995	996	997	998	999	1000

An example of the exported Sample list is shown as follows:

	A	B	C	D	E	F	G	H
1		X						
2		1 Accession						
3		2 Peptide Sequence						
4		3 FrameID						
5		4 II_171103_IAMVung_00_REP4_cdm122_HCD_OT_val45_HCD30_01_4ul						
6		5 II_171103_IAMVung_00A_cdm122_HCD_OT_val45_HCD30_01_4ul						
7		6 II_171103_IAMVung_00B_cdm122_HCD_OT_val45_HCD30_01_4ul						
8		7 II_171103_IAMVung_00C_cdm122_HCD_OT_val45_HCD30_01_4ul						
9		8 II_171103_IAMVung_01A_cdm122_HCD_OT_val45_HCD30_01_4ul						
10		9 II_171103_IAMVung_01B_cdm122_HCD_OT_val45_HCD30_01_4ul						
11		10 II_171103_IAMVung_01C_cdm122_HCD_OT_val45_HCD30_01_4ul						
12		11 II_171103_IAMVung_01D_cdm122_HCD_OT_val45_HCD30_01_4ul						
13		12 II_171103_IAMVung_02A_cdm122_HCD_OT_val45_HCD30_01_4ul						
14		13 II_171103_IAMVung_02B_cdm122_HCD_OT_val45_HCD30_01_4ul						
15		14 II_171103_IAMVung_02C_cdm122_HCD_OT_val45_HCD30_01_4ul						
16		15 II_171103_IAMVung_02D_cdm122_HCD_OT_val45_HCD30_01_4ul						
17		16 II_171103_IAMVung_03A_cdm122_HCD_OT_val45_HCD30_01_4ul						
18		17 II_171103_IAMVung_03B_cdm122_HCD_OT_val45_HCD30_01_4ul						
19		18 II_171103_IAMVung_03C_cdm122_HCD_OT_val45_HCD30_01_4ul						
20		19 II_171103_IAMVung_03D_cdm122_HCD_OT_val45_HCD30_01_4ul						
21		20 II_171103_IAMVung_04A_cdm122_HCD_OT_val45_HCD30_01_4ul						
22		21 II_171103_IAMVung_04B_cdm122_HCD_OT_val45_HCD30_01_4ul						
23		22 II_171103_IAMVung_04C_cdm122_HCD_OT_val45_HCD30_01_4ul						
24		23 II_171103_IAMVung_04D_cdm122_HCD_OT_val45_HCD30_01_4ul						
25		24 II_171103_IAMVung_05A_cdm122_HCD_OT_val45_HCD30_01_4ul						
26		25 II_171103_IAMVung_05B_cdm122_HCD_OT_val45_HCD30_01_4ul						
27		26 II_171103_IAMVung_05C_cdm122_HCD_OT_val45_HCD30_01_4ul						
28		27 II_171103_IAMVung_05D_cdm122_HCD_OT_val45_HCD30_01_4ul						
29		28 II_171103_IAMVung_06A_cdm122_HCD_OT_val45_HCD30_01_4ul						
30		29 II_171103_IAMVung_06B_cdm122_HCD_OT_val45_HCD30_01_4ul						
31		30 II_171103_IAMVung_06C_cdm122_HCD_OT_val45_HCD30_01_4ul						
32		31 II_171103_IAMVung_06D_cdm122_HCD_OT_val45_HCD30_01_4ul						
33		32 II_171103_IAMVung_07A_cdm122_HCD_OT_val45_HCD30_01_4ul						
34		33 II_171103_IAMVung_07B_cdm122_HCD_OT_val45_HCD30_01_4ul						
35		34 II_171103_IAMVung_07C_cdm122_HCD_OT_val45_HCD30_01_4ul						
36		35 II_171103_IAMVung_07D_cdm122_HCD_OT_val45_HCD30_01_4ul						
37		36 II_171103_IAMVung_08A_cdm122_HCD_OT_val45_HCD30_01_4ul						
38		37 II_171103_IAMVung_08B_cdm122_HCD_OT_val45_HCD30_01_4ul						
39		38 II_171103_IAMVung_08C_cdm122_HCD_OT_val45_HCD30_01_4ul						
40		39 II_171103_IAMVung_08D_cdm122_HCD_OT_val45_HCD30_01_4ul						
41		40 II_171103_IAMVung_09A_cdm122_HCD_OT_val45_HCD30_01_4ul						
42		41 II_171103_IAMVung_09B_cdm122_HCD_OT_val45_HCD30_01_4ul						
43		42 II_171103_IAMVung_09C_cdm122_HCD_OT_val45_HCD30_01_4ul						
44		43 II_171103_IAMVung_09D_cdm122_HCD_OT_val45_HCD30_01_4ul						
45		44 II_171103_IAMVung_10A_cdm122_HCD_OT_val45_HCD30_01_4ul						
46		45 II_171103_IAMVung_10B_cdm122_HCD_OT_val45_HCD30_01_4ul						
47		46 II_171103_IAMVung_10C_cdm122_HCD_OT_val45_HCD30_01_4ul						
48		47 II_171103_IAMVung_10D_cdm122_HCD_OT_val45_HCD30_01_4ul						

II. Protein quantification

1. Input files & load IonStarStat

Annotated frame list (.csv) A list of filtered frames with annotation (peptide sequence and protein AC).

Grouped sample list (.txt) A list containing all sample names and corresponding group information.

To generate the group list for quantification, open the sample list file (.csv) in excel. The file should look like this:

	A	B
1	Rawfiles	GroupID
2	II_171103_JAVIung_00_REP4_dtm122_HCD_OT_val45_HCD30_01_4ul	A
3	II_171103_JAVIung_00A_dtm122_HCD_OT_val45_HCD30_01_4ul	A
4	II_171103_JAVIung_00B_dtm122_HCD_OT_val45_HCD30_01_4ul	A
5	II_171103_JAVIung_00C_dtm122_HCD_OT_val45_HCD30_01_4ul	A
6	II_171103_JAVIung_01A_dtm122_HCD_OT_val45_HCD30_01_4ul	B
7	II_171103_JAVIung_01B_dtm122_HCD_OT_val45_HCD30_01_4ul	B
8	II_171103_JAVIung_01C_dtm122_HCD_OT_val45_HCD30_01_4ul	B
9	II_171103_JAVIung_01D_dtm122_HCD_OT_val45_HCD30_01_4ul	B
10	II_171103_JAVIung_02A_dtm122_HCD_OT_val45_HCD30_01_4ul	C
11	II_171103_JAVIung_02B_dtm122_HCD_OT_val45_HCD30_01_4ul	C
12	II_171103_JAVIung_02C_dtm122_HCD_OT_val45_HCD30_01_4ul	C
13	II_171103_JAVIung_02D_dtm122_HCD_OT_val45_HCD30_01_4ul	C
14	II_171103_JAVIung_03A_dtm122_HCD_OT_val45_HCD30_01_4ul	D
15	II_171103_JAVIung_03B_dtm122_HCD_OT_val45_HCD30_01_4ul	D
16	II_171103_JAVIung_03C_dtm122_HCD_OT_val45_HCD30_01_4ul	D
17	II_171103_JAVIung_03D_dtm122_HCD_OT_val45_HCD30_01_4ul	D
18	II_171103_JAVIung_04A_dtm122_HCD_OT_val45_HCD30_01_4ul	E
19	II_171103_JAVIung_04B_dtm122_HCD_OT_val45_HCD30_01_4ul	E
20	II_171103_JAVIung_04C_dtm122_HCD_OT_val45_HCD30_01_4ul	E
21	II_171103_JAVIung_04D_dtm122_HCD_OT_val45_HCD30_01_4ul	E
22	II_171103_JAVIung_05A_dtm122_HCD_OT_val45_HCD30_01_4ul	F
23	II_171103_JAVIung_05B_dtm122_HCD_OT_val45_HCD30_01_4ul	F
24	II_171103_JAVIung_05C_dtm122_HCD_OT_val45_HCD30_01_4ul	F
25	II_171103_JAVIung_05D_dtm122_HCD_OT_val45_HCD30_01_4ul	F
26	II_171103_JAVIung_06A_dtm122_HCD_OT_val45_HCD30_01_4ul	G
27	II_171103_JAVIung_06B_dtm122_HCD_OT_val45_HCD30_01_4ul	G
28	II_171103_JAVIung_06C_dtm122_HCD_OT_val45_HCD30_01_4ul	G
29	II_171103_JAVIung_06D_dtm122_HCD_OT_val45_HCD30_01_4ul	G
30	II_171103_JAVIung_07A_dtm122_HCD_OT_val45_HCD30_01_4ul	H
31	II_171103_JAVIung_07B_dtm122_HCD_OT_val45_HCD30_01_4ul	H
32	II_171103_JAVIung_07C_dtm122_HCD_OT_val45_HCD30_01_4ul	H
33	II_171103_JAVIung_07D_dtm122_HCD_OT_val45_HCD30_01_4ul	H
34	II_171103_JAVIung_08A_dtm122_HCD_OT_val45_HCD30_01_4ul	I
35	II_171103_JAVIung_08B_dtm122_HCD_OT_val45_HCD30_01_4ul	I
36	II_171103_JAVIung_08C_dtm122_HCD_OT_val45_HCD30_01_4ul	I
37	II_171103_JAVIung_08D_dtm122_HCD_OT_val45_HCD30_01_4ul	I
38	II_171103_JAVIung_09A_dtm122_HCD_OT_val45_HCD30_01_4ul	J
39	II_171103_JAVIung_09B_dtm122_HCD_OT_val45_HCD30_01_4ul	J
40	II_171103_JAVIung_09C_dtm122_HCD_OT_val45_HCD30_01_4ul	J
41	II_171103_JAVIung_09D_dtm122_HCD_OT_val45_HCD30_01_4ul	J
42	II_171103_JAVIung_10A_dtm122_HCD_OT_val45_HCD30_01_4ul	K
43	II_171103_JAVIung_10B_dtm122_HCD_OT_val45_HCD30_01_4ul	K
44	II_171103_JAVIung_10C_dtm122_HCD_OT_val45_HCD30_01_4ul	K
45	II_171103_JAVIung_10D_dtm122_HCD_OT_val45_HCD30_01_4ul	K

The column headers are **Rawfiles** and **GroupID** (case-sensitive), and samples in the same group should be assigned the same GroupID. The GroupID can be any combinations of alphabets and numbers.

Open **Run IonStar code.R** and run Line 1 to load IonStarStat.

2. Removal of redundant quantitative features (frames)

a. Open **Run IonStar code.R**.

b. Modify the following lines for input and output files:

Line 2~4

```
2 setwd("File directory")
3 rawfile <- "Raw_input_IonStar.csv"
4 condfile <- "Group list.txt"
```

Line 13~14

```
14 wrote.csv(quan,"Protein quantitative result.csv")
15 write.csv(exprs(cdata),"Peptide quantitative result.csv")
```

c. Run Line 5~7 to read the annotated frame list and the grouped sample list into R environment.

It is recommended to check the grouping information by entering **condition** in the R Console:

```
> condition
[1] A A A A B B B B C C C C D D D D E E E E F F F F G G G G H H H H I I I I J J J J K K K K
Levels: A B C D E F G H I J K
```

d. Run Line 8 to remove redundant frames (i.e. frames assigned to multiple peptide sequences) will also be removed to avoid ambiguity in quantification. The number of protein before and after removal, as well as the number of redundant frames removed will be reported in the console:

```
Input 6527 proteins.
19526 duplicated frames founded.
6465 proteins left after filtering.
```


3. *Frame-to-peptide aggregation and inter-sample normalization*

Run Line 9 to perform frame-to-peptide aggregation and inter-sample normalization of peptide intensities. Parameters include:

summarize = TRUE, FALSE (Whether to aggregate frame-level data to peptide level)

method = NULL, "TIC", "Quantiles" (Which normalization method is used, TIC being the default)

An example of sample intensities before and after normalization is shown below:



4. Peptide-level outlier detection

Run Line 10 to perform outlier peptide detection. IonStar uses Principal Component-based Outlier Detection (*PCOut*) for outlier detection, which is tailored for multi-condition comparison (≥ 3 conditions including control).

```
> cdata<-OutlierPeptideRM(ndata,condition,variance=0.7,critM1=1/3,critM2=1/4,ratio=TRUE)
13734 outliers were removed; 46268 peptides left after outlier removal.
```

Parameter **variation** = 0.7~0.9 can be adjusted according to the stringency needed for outlier detection. The higher the value the more outliers will be rejected.

For case-control comparison, please set parameter **ratio** = FALSE. Alternatively, Grubb's test can be used for outlier rejection. The Grubb's test function will be available in the next build of IonStarStat.

NOTE: during this step, errors could occasionally occur if there are too many missing values in certain frames. In this case, run Line 17~20 to remove frames with missing values, then restart the workflow from Line 8.

5. Shared peptide removal (Optional)

Run Line 11 to remove shared peptides (i.e. peptides inferred to ≥ 2 unique protein groups, a.k.a. degenerate peptides). This step is optional for protein quantification, as many highly abundant proteins share a large proportion of homologous sequence domains. Removal of these peptides could be counterproductive for quantification. However, in specific cases, such as quantification of mixed-species samples, removal of shared peptides with species ambiguity is necessary to obtain species-specific quantitative results.

6. Peptide-to-protein aggregation

Run Line 12 to obtain the quantitative results. Parameters include:

method = “sum”, “fit”

“sum” will simply add the intensities of all peptides inferred to the same protein to obtain protein quantitative intensities.

“fit” will apply General Linear Mixed Model (GLMM) to summarize peptide intensities to protein level.

7. Export quantitative results

Run line 13~14 to export peptide and protein quantitative results (.csv).

The output (Protein quantitative result.csv) contains raw quantitative intensities of all proteins as follows:

MODULE 4: POST-QUANTIFICATION DATA PROCESSING



StarGazer, the module for post-quantification data processing, will be made available in the next build of IonStar.

CONTACT INFORMATION

For questions, suggestions, and other topics about IonStar, feel free to contact us:

Shichen Shen shichens@buffalo.edu

Xue Wang xwang79@buffalo.edu

Jun Qu junqu@buffalo.edu